

# werner



Purified  
*Water*

Highly Purified  
*Water*

Concepts  
Development  
Engineering  
Construction  
Qualification  
Service

# Pharmaceutical *Water* . . . . .



## Quality and Diversity.

**The quality and diversity our technologies ensure detailed and total solutions to meet specific customer needs.**

This standard has enabled us to expand our business for pharmaceutical water systems in more than ten years to our strongest segment. As a result, Wilhelm Werner GmbH has become a leading provider for pharmaceutical water solutions, particularly for small and medium-sized companies.

Based on many years of experience in ultra clean processes, such as those in microelectronics and extensive engineering expertise, we are able to serve all fields of sterile production along with the conventional pharmaceutical industry, including food production, cosmetics manufacture, biotechnology and all other fields in the life science area. Production companies working upstream or parallel to these processes, such as suppliers of pharmaceutical packaging or research companies involved with scale-up, are also individually supplied by us with their own purification systems.



We understand how to harmonize customer expectations and requirements in an appropriate way: we always combine statutory requirements (such as EP, USP, JP) and FDA and GMP guidelines with individual needs for plant safety, space considerations and the latest technology. This has brought us great success so far, not the least


because we are a small, quality-focused, owner-operated company. Every customer always has a single point of contact in our team who knows what is needed. Flexibility, friendliness and special expertise are ordinary characteristics for all of our employees. Our customers appreciate this at the highest level.

## Total Customer Care.

All plant development processes including after-sales service are documented and certified: since 1996 a quality management system in accordance with DIN EN ISO 9001:2000 has been in successful use. Customer audits by major companies have continuously verified the effectiveness and reliability of the QM system.

We provide you with support for each phase of your project and guarantee the best implementation and comprehensive service.

- Expert advice
- Determination of basic requirements
  - Basic engineering
- Individual concept development
- Detailed engineering in accordance to FDA and cGMP guidelines
  - Individual creation of test documentation (DQ, IQ, OQ)
- Design qualification (DQ)
  - Plant construction
  - Automation
  - Installation and FAT
  - Operational startup
- Installation qualification (IQ)
- Operational qualification (OQ)
  - SAT and training
- Documentation for later validation
- Support during PQ
  - Customer service
- Maintenance service
  - Spare parts service
  - Hotline service
  - Calibration service
  - After-sales service
- Plant and process optimization
  - Upgrade and redesign of existing plants



# Purified Water . . . . .



## Purified Water Made to Measure. Custom Designed WERNER PW Systems.

Each plant for the production of purified water is custom work, whether its capacity is 100 l/h or over 25 m<sup>3</sup>/h. After all, each customer has individual system requirements: An SME engaging in contract manufacture has needs which differ from an international conglomerate. Thus no plant is entirely like any other. All WERNER PW basic systems meet EP, USP and JP guidelines and ensure the production process, product quality and profitability. The GMP-compliant design confirms to FDA and ISPE requirements.



## PW Basic System.

Ready-to-connect, pretested compact systems for purified water production. The compact, package unit design on a stainless steel frame offers all purification technologies, which are specifically adapted to the present tap water quality, such as

- automatic backwash pre-filters
- pipe separators
- duplex softening units, serially connected, quality-controlled, optional sanitisation with hot water > 80°C
- alternative conditioning
- UV+Plus recirculation system
- one- or two-pass reverse osmosis
- membrane degassing for CO<sub>2</sub> reduction
- CEDI continuous electrodisinfection
- chemical sanitisation or optional hot water sanitisation > 80°C
- pipework made of PVDF-HP by BCF welding technology or AISI 316L orbital welded
- state-of-the-art PLC technology with Siemens S7 PLC, touch panel or panel PC
- data storage in compliance with 21 CFR Part 11

## Custom Work in the Details. Custom Overall.

We are the right partner for large production capacities even well above 25 m<sup>3</sup>/h as well as for difficult conditions, such as intermittent use, seasonal variations in production, reconstruction projects with ongoing production or critical tap water. All the expertise gained from customer-specific solutions and the latest developments in technology contribute to creating the individual WERNER PW system solution for every requirement.

- Pre-treatment**
  - alternative preconditioning technologies
  - serial softening systems with chemical or hot water sanitisation
  
- Recirculation system UV+Plus**
  - continuous circulation of RO and CEDI
  - 100 % recycling of CEDI wastewater in all operating modes
  - 100 % recycling of RO concentrates in the recirculation mode
  - UV disinfection in the RO feed (1,250 to 3,000 J/m<sup>3</sup>)
  
- Sanitisation**
  - dialog-supported chemical sanitisation program
  - optional hot water sanitisation > 80 °C
  
- Expandable plants**
  - modular system concept
  - subsequent capacity expansion up to 100 %
  
- EEx**
  - systems for EEx Zone 1 with EEx-p encapsulation



# Highly Purified Water . . . . .



## Highly Purified Water. Costs to Fit Your Budget.

The inclusion of the quality “highly purified water” (HPW) in the European Pharmacopeia on June 1, 2002 enables a new plant design with an additional depyrogenation step after the purified water production. This allows considerable savings of capital expenditures and operating costs for the final rinse of parenteral containers and equipment which previously had to be cleaned with WFI.



## HPW Basic System.

Ready-to-connect, pretested compact systems for highly purified water or WFI (USP/JP) production. The compact, package unit design on a stainless steel frame offers all purification technologies, which are specifically adapted to the tap water quality present, such as

- automatic backwash pre-filters
- pipe separators
- duplex softening units, serially connected, quality-controlled, optional sanitisation with hot water > 80°C
- alternative conditioning
- UV+Plus recirculation system
- one- or two-pass reverse osmosis
- membrane degassing for CO<sub>2</sub> reduction
- CEDI continuous electrodisinfection
- ultrafiltration, 6,000 Dalton cut-off
- chemical sanitisation or optional hot water sanitisation > 80°C
- pipework made of PVDF-HP by BCF welding technology or AISI 316L orbital welded
- state-of-the-art PLC technology with Siemens S7 PLC, touch panel or panel PC
- data storage in compliance with 21 CFR Part 11



## Custom Work in the Details. Custom Overall.

Werner HPW systems can be combined with existing PW production systems or will be designed directly as a complete HPW system. The unit includes all the important purification steps for producing HPW; entire design and construction are in accordance with all guidelines and recommendations of the FDA, GAMP and ISPE. Based on the Werner PW basic systems, Werner provides HPW systems with the following features:

- sealless ultrafiltration unit with a cut-off of 6,000 Dalton
- optional recirculation pump
- optional hot water sanitisation at  $> 80^{\circ}\text{C}$
- compact, modular construction, also as a single-stage unit downstream from existing PW systems
- optional EEx execution



# Pharmaceutical Water.....



## Storage.

After the purification of PW or HPW, engineering of the customer-specific storage and distribution system takes place so product water can be delivered to each point of use in the sterile distribution system without recontamination. Here too we consider the customer's specific conditions and produce an individually engineered design.

- 316L SS storage tanks for PW and HPW
- frequency controlled distribution pump
- optional dual pump station, continuous flow-through
- UV disinfection / ozone reduction
- tap management
- uncompensated conductivity measurement
- flow management
- ozonation
- online TOC measurement
- optional thermal sanitisation at > 80 °C
- cooling with DTS heat exchanger



## Distribution.

PVDF-HP has established itself as a material for all distribution systems particularly because of its high surface quality, with an average  $Ra < 0.20 \mu\text{m}$ , and low ion and TOC leach-out. With its temperature resistance up to  $140^\circ\text{C}$ , a corresponding piping calculation can be performed and permanent hot storage at  $85^\circ\text{C}$  or periodic heat sanitisation implemented.


Likewise the entire distribution system can be made of 316L stainless steel, with the option to improve the surface quality from  $Ra < 0.80 \mu\text{m}$  to  $Ra < 0.40 \mu\text{m}$ .



## Automation.

Planning and producing of the entire automation system requires experience and knowledge of the overall process. A team consisting of a project leader and automation engineer supports the project from the design phase through its conclusion, enabling this sensitive field to be covered with complete competence.

- project coordination
- one-stop implementation
- hardware design
- software design
- control technology system design
- software validation
- design and implementation with EPLAN/WSCAD, ProTool/WinCC







# Pharmaceutical Water.....

## Validation is Teamwork.

Together with the customer, we develop a basis for the entire qualification. Each customer has different benchmarks for the internal validation process.

We accept the requirements of the qualification management, couple these with our experiences and in this manner define a common path for the validation process right at the start of the project based on teamwork with the customer from the beginning.

The validation of the PW and HPW systems is performed in accordance with FDA, cGMP, GAMP, USP and PharmEUR regulations as well as ISPE baselines.



## Service is Teamwork.

As we understand it, customer relations even after operational startup is also teamwork, because only thorough instruction and training of operating personnel ensure high customer satisfaction and a long service life for the system.

The associated maintenance service is defined in parallel to the time for startup. Numerous optional services are offered as part of the standard basic and primary maintenance; these can also be combined in a complete service package:

- basic and primary maintenance as a customer-specific service
- user training of operating personnel
- critical tap water analysis by in-house or contract laboratory
- calibration service for all measurement points subject to calibration
- certified reference measurement using ASTM standards
- provision and management of system-specific test equipment
- remote maintenance via modem
- development of operation models
- customer-specific spare parts provision



## References

### international

Alkaloid AD · Skopje Macedonia  
 Arabio · Jeddah · Saudi Arabia  
 AC Helcor · Baia Mare Romania  
 Alsifcom · Cluj-Napoca Romania  
 Ardea Pharma · Sevetin Czechoslovakia  
 Beiersdorf Thailand Co. Ltd. · Bangkok Thailand  
 Esco France S.A. · Dombasle France  
 Farmex Companx S.R.L. · Bucharest Romania  
 Pharmex Rom Industry S.R.L. · Ilfov Romania  
 Farmacon S.R.L. · Brasov Romania  
 FDS Pharma ASS · St. Petersburg Russia  
 Galax d.o.o. · Murska Sobota Slovenia  
 Gematek OOO · Tver Russia  
 Hoechst Biotika spol.sr.o. · Martin Slovakia  
 ICN-Farm AD · Skopje Macedonia  
 ICN Czech Republic A.S. · Roztoky Czechoslovakia  
 Institute of Virology · Vaccines and Sera Torlak · Belgrade Serbia  
 Institute Apicol · Bucharest Romania  
 Medipharm (Pvt) Ltd · Lahore Pakistan  
 MTH (Binnopharm) · Zelenograd Russia  
 Nivea (Shanghai) Company Ltd. · Qing Pu China  
 PharmaConstruct Ltd. · Orgsintez Russia  
 Shifa Pharmaceutical Industries · Aleppo Syria  
 University Pharmacia · Prag Czechoslovakia

### national

Altana Pharma AG · Konstanz, Singen  
 Ascalon GmbH · Berggießhübel  
 BASF AG · Minden  
 Baxter Oncology GmbH · Halle Westf.  
 Bayer AG · Dormagen  
 Bayer Bitterfeld GmbH · Greppin  
 Bela-pharm GmbH & Co. KG · Vechta  
 Boehringer Ingelheim Pharma GmbH & Co. KG · Biberach  
 Boehringer Ingelheim MicroParts GmbH · Dortmund  
 Braun GmbH · Bahlingen  
 C.P.M. Contract Pharma GmbH · Feldkirchen  
 Chemische Fabrik Lehrte – Dr. Andreas Kossel GmbH · Lehrte  
 Dade Behring Marburg GmbH · Marburg  
 Dologiet Arzneimittel GmbH & Co. KG · St. Augustin  
 Dr. August Wolff GmbH · Bielefeld  
 Dr. R. Pfleger GmbH · Bamberg  
 Excella GmbH · Feucht  
 GeneScan AG · Freiburg  
 Girindus AG · Halle  
 Hexal AG · Holzkirchen  
 Linde AG, Dresden  
 MAQUET Cardiopulmonary AG · Hechingen  
 Medtronic Sfamor Daneck · Deggendorf  
 Membrana GmbH · Obernburg  
 Merck KgaA · Darmstadt  
 Mibe GmbH · Brehna  
 NextPharma, Pharbil Waltrop GmbH · Waltrop  
 Nordmark Arzneimittel GmbH & Co. KG · Uetersen  
 Norbítec GmbH · Uetersen  
 Organon Development GmbH · Waltrop  
 Oxoid Biotech GmbH · Wesel  
 Pegasus Pharma GmbH · Hannover  
 Pharmazeutische Fabrik Lichtenheldt GmbH · Wahlstedt  
 Qiagen GmbH · Hilden  
 Roche Diagnostics GmbH · Mannheim  
 Roha Arzneimittel GmbH · Bremen  
 Rottendorf Pharma GmbH · Enningerloh  
 Sandoz Pharmaceuticals GmbH · Gerlingen  
 Sanorell-Pharma GmbH · Bühl  
 Solvay Pharmaceuticals GmbH · Nienburg  
 Temmler Pharma GmbH & Co. KG · Marburg  
 Trion Pharma GmbH · München  
 Tutogen Medical GmbH · Neunkirchen  
 Wacker Chemie GmbH · Nünchritz  
 Wala Heilmittel GmbH · Bad Boll



## Uncompromising in Quality. Consequent in Details.

### **Experience and innovation.**

For more than two decades our name has stood for the highest quality in the purification of pure and ultrapure water. We are in demand as a partner internationally with dedicated sales organizations and subsidiaries in many European countries. With a focus on quality, an understanding of market needs, innovative drive and flexibility we have grown to become one of the leading providers of pure and ultrapure water purification systems on the German and international markets.

### **Products and service from a single source.**

Increasingly precise and challenging analytical and production technologies in research and industry require innovative water purification methods.

We offer everything – from the smallest laboratory system with 0.5 l/min up to pilot plant systems with 25 m<sup>3</sup>/h: precise, customized solutions from a single supplier. This is backed up by a highly motivated team which develops and produces practical solutions and covers the range from analysis and consulting to CAD design, in-plant manufacture, operational startup and maintenance.

### **Certified quality.**

Since 1996 we have been certified as meeting the highest internationally recognized quality requirements. WERNER was one of the first providers of water purification systems to achieve ISO certification. This established a basis for future security for companies, customers and business partners.

Subject to technical changes without notice.

Last revised 04/2009

# werner

Wilhelm Werner GmbH  
Reinstwassertechnik  
Maybachstraße 29  
D-51381 Leverkusen  
tel +49 2171 7675-0  
fax +49 2171 7675-10  
e-mail: [info@werner-gmbh.com](mailto:info@werner-gmbh.com)  
[www.werner-gmbh.com](http://www.werner-gmbh.com)